


# CO<sub>2</sub> Laser

C-Las



# CO<sub>2</sub> Laser

## › Laser cutting



### Flexible and Versatile

Maximum use is assured with carefully designed accessories for many applications.

**Improved surgical results:** The tactile sensation of the free style hand piece combined with the small spot size enable excellent treatment results with minimal recovery time.

## › Skin ablation



### Sure Scan Ablation

Increased precision and adjustable exposure time make this the laser of choice for alisar la piel, Wrinkle Removal, Pigmented Lesions, etc.

## › C-Las with articulated arm



### C-Las also offers an optional

Articulated arm delivery system for Colposcopy and Micromanipulators.

## Technical data

Laser	CO <sub>2</sub> , sealed ceramic tube
Wavelength	10600 nm
Output power	Min. 30 Watt at laser head
Beam profile	Quasi Gauss at end of Wave Guide
Operating mode	CW, and pulsed
Beam guiding	Flexible hollow wave guide 1.8 mm OD or Articulated arm
Length of Wave Guide	~ 2.6 m, different length upon request
Accessories	Universal hand piece, Focus hand piece, Flying spot scanner, Micromanipulator
Dimensions	478 x 200 x 491 mm
Weight	18 kg
Mains supply	100-260 V, 50/60 Hz, 10 A

Datas are subject to change, to maintain the highest standard



**A.R.C. Laser**  
certifies that the product complies to  
21 CFR 1040.10 and 1040.11  
EN 13485/2003 - 93/42 EWG