

# Slit Lamp

PCL3



# Excellent Price



The slit lamp is an instrument which is routinely used in every office of an Ophthalmologist. Quality and easy use combined with good optics is the advantage of our PCL3.

All attachments -like video -are available.

In combination with our **Classic KTP Laser**, this is the instrument of choice for Photocoagulation.

## Feather Touch Joystick

Longitudinal (in/out)	90 mm
Lateral (left/right)	100 mm
Vertical (up/down)	30 mm
Chin Rest Range	80 mm
Voltage	Universal AC input 110V-240V (50/60 Hz)



## Bright Light Illumination

Slit Length	0-14 mm
Slit Apertures	0.3 / 5.5 / 9 / 14 mm
Slit Rotation	0-180°
Filters	Red free, cobalt blue, heat absorbing
Illumination	6V 20W halogen



## Video Attachment Available

Model	PCL3
Microscope	Galilean
Magnification	3 step rotation; 10x / 16x / 25x
Eyepiece	Wide field 12.5x
PD Range	48.5-80 mm, Adjustment $\pm 6$ Diopter