

Accessories

Version: 01-2008

Page: 1 / 4

KTP Laser 532nm High Power 8 Watt NUVOLAS

for fiber Ø 200µm
 Ø 300µm
 Ø 600µm











- Vascular lesions
- Pigmented lesions

Accessories

Version: 01-2008

Page: 2 / 4






Item-No.	Recommended Accessories Dermatology	Picture
LL13002	Fiber orange, 2 click connectors Length: 1,65m (for focus or bleaching handpiece) 300µm fiber (connector white)	
HS11008 HS11009 HS11010	Focus handpiece (for 300µm fiber LL13002) Ø 0,25mm (blue) Ø 0,60mm (red) Ø 1,20mm (green)	





Item-No.	Recommended Accessories Surgery	Picture
LL13003	Barefiber orange, disposable 1 click connector, Length: 1,65m 200µm fiber (connector red)	
LL13001	Barefiber orange, disposable 1 click connector, Length: 1,65m 300µm fiber (connector white)	
LL13008	Barefiber orange, sterilisable 1 click connector, Length: 1,65m 600µm fiber (connector green)	
HS11013	Surgery handpiece (in combination with a barefiber)	
ZU01024	Stainless steel cannula (in combination with surgery handpiece) curved 2cm, outer Ø 1,0mm	
ZU01016	Stainless steel cannula (in combination with surgery handpiece) curved 5cm, outer Ø 1,0mm	

Accessories

Version: 01-2008

Page: 3 / 4


Item-No.	Recommended Accessories Surgery	Picture
ZU01017	Stainless steel cannula (in combination with surgery handpiece) straight 5cm, outer Ø 1,0mm	
ZU01031	Stainless steel cannula (in combination with surgery handpiece) straight 5cm, outer Ø 0,7mm	
ZU01030	Stainless steel cannula (in combination with surgery handpiece) straight 8cm, outer Ø 0,7mm	
ZU01027	Plastic tip curved, disposable, unit pack: 10 pieces	
WE01275	Sterilisation box Size: 152x51x19mm	

Item-No.	Further Accessories	Picture
WE01209	Stripper, adjustable	
WE01003	Cutter	
AS01003	Safety goggle 488-532nm	 <small>Pic. similar</small>
AS01016	Safety goggle for patient Spectra Shield, black	 <small>Pic. similar</small>

Accessories

Version: 01-2008

Page: 4 / 4

Item-No.	Further Accessories	Picture
AS01012	Eye shield, size: middle Stainless steel, surface: matt	 Pic. similar
KB01137	Connecting cable: from laser to remote control	