


# VARIO CLASSIC

## ND:YAG + ARGON-LASER

**IDEAL IN  
CLINIC &  
OFFICE**









**Capsulotomy**  
**Iridotomy**  
**Retina**  
**Coagulation**






### VARIO Q-LAS + CLASSIC 514

Amount	Order #	Description
1 x	LS09002	Q-LAS: Q-switched Nd:YAG Laser 1064 nm, <sup>2</sup> slit lamp PCL5, internal eye safety filter, parallel binocular tube
	LS04100	CLASSIC 514, Diode pumped laser 514 nm, <sup>1</sup> fiber 62 μm SMA plug Slit lamp PCL5, convergent binocular tube, Super View Adapter
	SL11050	VARIO table (116 cm width) <sup>3</sup>
		2 x dust protection cover


### ACCESSORIES FOR MORE COMFORT

<p>PCS=1   ZU11038</p>  <p>Elbow rest, 24 mm</p>	<p>PCS=1   ZU11039</p>  <p>Elbow rest, 48 mm</p>
<p>PCS=1   WE01243</p>  <p>Iridotomy lens, Abraham</p>	<p>PCS=1   WE01244</p>  <p>Capsulotomy lens, Abraham</p>
<p>PCS=1   WE01242</p>  <p>3 mirror universal Gonio lens, 18 mm</p>	<p>PCS=1   WE01245</p>  <p>NMR Mainster (standard) focal/grid</p>

### STANDARD DEVICE FEATURES

 Hand trigger	 Fixed stand with sliders	 LED illumination
---	---	---



### EXCLUSIVE UPGRADE FOR Q-LAS (EXTRA CHARGE)

<p>PCS=1   BG09093</p>  <p>Head-up Display</p>	<p><b>displays:</b></p> <ul style="list-style-type: none"> <li>• selected energy per pulse</li> <li>• total number of exposures</li> </ul>
--	--

### DEVICE OPTION - ON REQUEST (EXTRA CHARGE)

 Stand with castors (PCS=4)	 BG04008 Foot switch
---	--

### EYE SAFETY

	<p>AS01033 safety goggle PCS=1</p> <p>445/810/980 1064 nm</p>
	<p>AS01003 safety goggle PCS=1</p> <p>445/514/532/ 980/1064 nm</p>