


VARIO SLT

ND:YAG + SLT-LASER CITO

**IDEAL IN
CLINIC &
OFFICE**








**Capsulotomy
Iridotomy
Glaucoma
treatment**





Q-LAS

Amount	Order #	Description
1 x	LS09002	Q-LAS: 2 Q-switched Nd:YAG Laser 1064 nm, slit lamp PCL5, internal eye safety filter, parallel binocular tube
	LS09010	CITO 532: 1 Frequency doubled, Q-switched Laser 532 nm, up to 10Hz of repetition rate, internal eye safety filter, slit lamp PCL5, parallel binocular tube
	SL11050	VARIO table (116cm width) 3 2 x dust protection cover


ACCESSORIES FOR MORE COMFORT

PCS=1 ZU11038  Elbow rest, 24 mm	PCS=1 ZU11039  Elbow rest, 48 mm
PCS=1 WE01243  Iridotomy lens, Abraham	PCS=1 WE01244  Capsulotomy lens, Abraham
PCS=1 WE01356  SLT Gonio Lens, Latina	

STANDARD DEVICE FEATURES

 Hand trigger	 Fixed stand with sliders	 LED illumination	 Parallel binocular tube
--	---	---	--



EXCLUSIVE UPGRADE (EXTRA CHARGE)

PCS=1 BG09093  Head-up Display	displays: <ul style="list-style-type: none"> • selected energy per pulse • total number of exposures
--	---

DEVICE OPTION - ON REQUEST (EXTRA CHARGE)

 Stand with castors (PCS=4)	 BG04008 Foot switch
---	--

EYE SAFETY

AS01033  safety goggle 445/810/980 1064 nm PCS=1
AS01003  safety goggle 445/514/532/ 980/1064 nm PCS=1