


**A.R.C.
LASER**

CLASSIC 514 + L-ord PHOTOCOAGULATION

**IDEAL IN
CLINIC &
OFFICE**



Ophthalmology
Retina
Coagulation

CLASSIC 514

Amount	Order #	Description
1 x	LS04100	CLASSIC 514 ① Diode pumped Laser 514 nm
	ZU12200	L-ord ② (retractable zoom-optic) integrated eye safety filter 514nm

ACCESSORIES FOR MORE COMFORT

PCS=1 ZU11038 Elbow rest, 24 mm	PCS=1 ZU11039 Elbow rest, 48 mm
PCS=1 WE01242 3 mirror universal Gonio lens, 18 mm	PCS=1 WE01245 NMR Mainster (standard) focal/grid
PCS=5 ME06895 Hygienic shield	

STANDARD DEVICE FEATURES

Holder for L-ord ③	BG06576 Adapter for tonometer or L-ord	BG04008 Foot switch	Transportation cases for L-ord and CLASSIC

EYE SAFETY

	AS01003 safety goggle	PCS=1
	445/514/532/ 980/1064nm	