


**A.R.C.
LASER**

COBRA

Nd:YAG + SLT-LASER

**IDEAL IN
CLINIC &
OFFICE**



**Capsulotomy
Iridotomy
Glaucoma
treatment**



VARIO SLT		
Amount	Order #	Description
1 x	LSXX	COBRA: Q-switched Nd:YAG Laser 1064 nm, Q-LAS slit lamp PCL5, internal eye safety filter, parallel binocular tube + Frequency doubled, Q-switched Laser 532 nm, up to 10Hz of repetition rate, internal eye safety filter, slit lamp PCL5, parallel binocular tube
	SL11040	Single table (86 cm width) ②
		2 x dust protection cover

ACCESSORIES FOR MORE COMFORT	
PCS=1 ZU11038 Elbow rest, 24 mm	PCS=1 ZU11039 Elbow rest, 48 mm
PCS=1 WE01243 Iridotomy lens, Abraham	PCS=1 WE01244 Capsulotomy lens, Abraham
PCS=1 WE01356 SLT Gonio Lens, Latina	PCS=5 WE08956 Hygienic shield

STANDARD DEVICE FEATURES			
Hand trigger	Fixed stand with sliders	LED illumination	Parallel binocular tube

EYE SAFETY		
	AS01033 safety goggle 445/810/980 1064 nm	PCS=1
	AS01003 safety goggle 445/514/532/ 980/1064 nm	PCS=1

DEVICE OPTION - ON REQUEST (EXTRA CHARGE)	
PCS=4 ME06340 Double castors	BG04008 Foot switch